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- Bernard A, Boumsell L, Dausset J, Milstein C and SchlossmanSF,eds:1984. Leukocyte Typing. New York: Springer-Verlag.p.28,41-42,44,196.
- McMichael AJ, ed:1987. Leukocyte a Typing III. Oxford: Oxford University Press.p.38,40,42,43,116,167,170-172,176,199,202,206,302-308,315,475.
- Reinherz EL, Haynes BF, Nadler LM and Bernstein ID:1986, Leukocyte Typing II. New York: Springer-Verlag. Vol 2, p.8, 15-20, 37.
- Reinherz EL and Schlossman SF:1980. The differentiation and function of human T lymphocytes. Cell 19:821-827.
- Aiuta F, Cerottini J-C, Coombs RRA, Cooper M, Dickler HB, Froland S, Fudenberg HH, Greaves MF, Grey HM, Kunkel HG, Natvig J, Preud'homme J-L, Rabellino E, Ritts RE, Rowe DS, Seligmann M,Siegal FP, Stjernsward J, Terry WD and Wybran J:1975. Identification, enumeration and isolation of B and T lymphocytes from human peripheral blood. Intemational Union of Immunological Scienties (IUIS), Report-July 1974. Clin Immunol and Immunopathol 3:584-597.
- Foon KA and Todd RF:1986. Immunologic classification of leukemia and lymphoma. Blood 68:1-31 .
- Drexler HG, Gignac SM and Minowada J:1988. Routine immunophenotyping of acute leukemias. Blut 57:327-339.
- Caligiuri M, Murray C, Buchwald D, Levine H,Cheney P, Peterson D, Komaroff AL and Ritz J:1987. Pheno-typic and functional deficiency of natural killer cells in patients with chronic fatigue syndrome. J Immnol 139:3306-3313.
- Reinherz EL, Meuer S, Fitzgerald KA, Hussey RE, Levine H and Schlossman SF:1982. Antigen recong- by human T lymphocytes is linked to surface expression of the T3 molecular complex. Cell 30: 735-743.

10. Meuer SC, Acuto O, Hussey RE, Hodgdon JC, Fitzgerald KA, Schlossman SF and Reinherz EL:1983. Evidence for the T3-associated 90K heterodimer as the T-cell antigen receptor. *Nature* 303:808-810.
11. Nadler LM, Anderson KC, Marti G, Bates M, Park E, Daley JF and Schlossman SF: 1983. B4, a human B lymphocyte-associated antigen expressed on normal, mitogen-activated, and malignant B lymphocytes. *J. Immunol.* 131:244-250.
12. Benjamini E and Leskowitz S: 1991. Disorders of the immune response. In: *Immunology: A Short Course. Second Edition.* New York: Wiley-Liss. p.211-244.
13. Reinherz EL, O'Brien C, Rosenthal P and Schlossman SF: 1980. The cellular basis for viral-induced Immunodeficiency: Analysis by monoclonal antibodies. *J Immunol* 125: 1269-1274.
14. Felsenstein D, Carney WP, Iacoviello VR and Hirsch MS: 1985. Phenotypic properties of atypical lymphocytes in cytomegalovirus-induced mononucleosis. *J Infect Dis* 152: 198-203.
15. Rinaldo CR, Ho M, Hamoudi WH, Gui X and DeBiasio RL: 1983. Lymphocyte subsets and natural killer cell responses during cytomegalovirus mononucleosis. *Infect Immun* 40: 472-477.
16. Goldstein G, Lifter J and Mittler R: 1982. Immunoregulatory changes in human disease detected by monoclonal antibodies to T lymphocytes. In: *Monoclonal Antibodies in Clinical Medicine.* McMichael AJ and Febre JW, eds. New York, NY: 1 Academic Press. p.39-70.
17. Schmidt RE: 1989. Monoclonal antibodies for diagnosis of immunodeficiencies. *Blut* 59: 200-206.
18. de Martini RM and Parker JW: 1989. Immunologic alterations in human immunodeficiency virus infection: A review. *J Clin Lab Anal* 3: 56-70.
19. Fauci AS: 1988. The human deficiency virus: Infectivity and mechanism of pathogenesis. *Science* 239: 617-622.
20. Taylor MGJ, Fahey JL, Detels R and Giorgi JV: 1989. CD4 percentage, CD4 number and C04:CD8 ratio in HIV infection: How to choose and how to use. *J AIDS* 2: 114-124.
21. Pedrazzini A, Freedman AS, Andersen J, Heflin L, Anderson K, Takvorian T, Canellos GP, Whitman J, Coral F, Ritz J and Nadler LM: 1989. Anti-B cell monoclonal antibody purged autologous bone marrow transplantation for B-cell non-Hodgkin's lymphoma: Phenotypic reconstitution and B-cell function. *Blood* 74: 2203-2211.
22. Preijers FWMB, DeWitte T, Wessels JMC, DeGast GC, Van Leeuwen E, Capel PJA and Haanen C: 1989. Autologous transplantation of bone marrow purged in vitro with anti-CD7-(WT1-) Ricin A Immunotoxin in T-cell lymphoblastic leukemia and lymphoma. *Blood* 74: 1152-1158.
23. Ramos EL, Turks LA, Leggat JE, Wood IG, Milford EL and Carpenter CB: 1989. Decrease in phenotypically defined T helper inducer cells (T4+4B4+) and increase in T suppressor effector cells (T8+2H4+) in stable renal allograft recipients. *Transplantation* 47: 465-471.
24. Beverly PCL and Callard RE: 1981. Distinctive functional characteristics of human 'T' lymphocytes defined by E rosetting or a monoclonal anti-T cells antibody. *EurJ Immunol* 11: 329-334.
25. Loken MR, Brosnan JM, Bach BA and Ault KA: 1990. Establishing optimal lymphocyte gates for immunophenotyping by flow cytometry. *Cytometry* 11 : 453-459.
26. Gebel HM, Lebeck LL, Jensik SC, Landay AL and Bray RA: 1989. Discordant expression of CD3 and T-cell receptor antigens on lymphocytes from patients treated with OKT3. *Transplantation Proceedings* 1: 1745-1746.
27. Tunnacliffe A, Olsson C, Traunacker A, Krissensen GW, Karjalainen K and De Le Heta A: 1989. The majority of CD3 epitopes are conferred by the epsilon chain. In *Leukocyte Typing IV, White Cell Differentiation Antigens.* W Knapp et al. eds, Oxford: Oxford University Press. p. 295-296.
28. Schroeder TI and Chatenoud L: Immunological monitoring during treatment with OKT3. Presented under the auspices of the American Society of Transplant Surgeons. Citation to be filled in when available from authors.
29. Koepke JA and Landay AL: 1989. precision and accuracy of absolute lymphocyte counts. *Clin Immunol Immunopathol* 52: 19-27.

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